

SMART



**SPIROBANK™
SMART**



Spirobank™ Smart



iSpirometry

Ultra portable, Real Time test


Suitable from 5 to 93 years!

Spirobank Smart is a complete inexpensive spirometer ideal for screening

Share results in pdf (WhatsApp, email etc.) and print directly



Patient Incentive



On-screen Results



Spirometry test:

PEF, FVC, FEV1, FEV1/FVC Ratio, FEV6 FEF 25/75



Available with both **DISPOSABLE** or **REUSABLE** single patient turbine flowmeter

Unparalleled Accuracy Unmatched Quality!

ISO 26782 Std (for spirometry parameters)

ISO 23747 Std (only for the Peak Flow parameter)

ATS/ERS 2019 standards