

Hypodermic Kit



DISPO VAN[®]

Sterile Single Use Needles

Manufactured using premium-quality cold-rolled stainless steel strips with the latest Japanese technology, Dispovan needles are ground on high-precision machines: assembled and packed on totally automated machines. Vision cameras ensure 100% automated quality control, enabling superior sharpness to minimize pain.

COMPONENTS:

NEEDLE HUBS : made of non-toxic, medical-grade Polypropylene compatible with the injectable medication.

CANNULAS : BURR FREE - 3 FACET LANCET POINT are made from SUS 304 stainless steel.

PACKING : The needles are packed in tamper-evident blister packs having easily peelable, medical-grade paper which ensure reliability of seal and product sterility.

PRECAUTIONS : Before use, check the integrity of the sealing of the pack. Store in clean Dry & insect free place. Avoid exposure to excessive heat. Avoid exposure to acidic/corrosive fumes, to avoid corrosion/rust.

Avoid piercing the needle through the wall in hard plastic container, it may result in needle blockage.

COLOUR CODE : As per ISO:6009.

SINGLE USE HYPODERMIC NEEDLE

COLOUR CODE AND SPECIFICATION

SIZE		COLOUR		LENGTH						
OD (mm)	Guage	Colour & Code		3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
0.40	27G	Grey	Grey							
0.45	26G	Brown	Brown							
0.50	25G	Orange	Orange							
0.55	24G	Medium purple	Medium purple							
0.60	23G	Blue	Blue							
0.70	22G	Black	Black							
0.80	21G	Green	Green							
0.90	20G	Yellow	Yellow							
1.10	19G	Creamish	Creamish							
1.20	18G	Pink	Pink							
1.60	16G	White	White							

• Conform to : ISO:7864/IS:10654 standards.

▲ Before use check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.

