

DISPO VAN[®]

Sterile Single Use Syringes

"Ribbon" Packed Dispovan, single use syringes without needle & needle mounted on top, are synonymous with quality

COMPONENTS:

BARRELS & PLUNGERS : made of non-toxic, medical-grade polypropylene compatible with the medication for short term contact drug delivery application.

The nozzle of the barrel is luer mount and 6% luer taper as per ISO:594.

GASKETS : Chemically inert and compatible with the medications for short term contact drug delivery application.

PACKING : DISPOVAN, single use syringes are individually ribbon packed in a double laminated plastic film. Ribbon packaging is superior to blister packaging due to higher puncture and water resistance, enabling extended shelf-life under adverse conditions. But for better container loadability we are also making transparent thermoformed peel able blister pack



PRODUCT RANGE		
Capacity	Type of Nozzle	Standard Needle Dimensions
1 ml.* (T.B.)	Luer mount with central nozzle	26G x 1/2" *0.45 x 13
2 ml.**	Luer mount/ central nozzle with	24G x 1", 23G x 1", 23G x 1 1/2", 22G x 1" *0.55 x 25 0.60 x 25 0.60 x 38 0.70 x 25 22G x 1 1/4", 22G x 1 1/2", 21G x 1", 21G x 1 1/2" *0.70 x 32 0.70 x 38 0.80 x 25 0.80 x 38
3 ml.	Luer mount/ central nozzle with	Specifications same as 2 ml
5 ml.**	Luer mount/ central nozzle with	Specifications same as 2 ml
10 ml.**	Luer mount with eccentric nozzle	22G x 1 1/4", 22G x 1 1/2", 21G x 1 1/2", 20G x 1 1/2" *0.70 x 32 0.70 x 38 0.80 x 38 0.90 x 38 Same as above
20 ml.	Luer mount with eccentric nozzle	Specifications same as 10 ml
50 ml.	Luer mount with eccentric nozzle	Without needle



**also available in Blister pack.

* Conform to : ISO:7886-1/ISO:10258 standards.

*All units are in mm

▲ Before use check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.



www.hmdhealthcare.com