



Improved
sharpness
for added
precision



GLASS VAN®

SCALPEL & SURGICAL BLADES

GLASS VAN®


SCALPEL & SURGICAL BLADES

GLASS VAN Surgical Blades

- Choice of Surgeons, Dentists, Chiropodists
- Made from premium quality Swiss High Carbon Steel
- Manufactured as per ISO : 7740
- CE Marked under Medical Device Directive (93/42/EEC)
- Individually packed in an easy peel aluminium pack
- Additionally protected for corrosion by VCI coated paper
- Sterilized by Gamma Radiation of 2.5 Mrad to ensure sterility assurance level of 10^{-6}

Product Range

Category	Commercial Designation No.
For small-size fittings (Handle No. 3)	10, 11, 12, 15
For large-size fittings (Handle No. 4)	20, 21, 22, 23, 24, 25



Packing Dimensions

Size of Shipping Carton (mm)	Primary Packaging (Pcs.)	Shipping Carton (Pcs.)	Gross Weight (Kg.)
540X180X180	100	5000	9.00
565X170X170	50	5000	9.00

GLASS VAN Surgical Scalpel with Plastic Handle

- Disposable Scalpel with high quality Stainless Steel Blades
- Ergonomically designed handle for better feel & precision
- Ethylene Oxide (EO) Gas sterilized
- Protected edge with easily removable plastic sheath
- Easy Peel Blister Pack in Medical Grade Paper complying with EN 868

Packing Dimensions

Size of Shipping Carton (mm)	Primary Packaging (Pcs.)	Shipping Carton (Pcs.)	Gross Weight (Kg.)
545X225X470	10	1000	12.20

