



Littmann®
Stethoscopes



L is for learning.

Quality makes a difference when you're learning the skill of auscultation to help you assess and monitor patients. Start with a 3M™ Littmann® Stethoscope for professional level acoustics. It's built to last and trusted by professionals worldwide.

The high acoustic sensitivity of 3M™ Littmann® Stethoscopes can help you learn to assess and monitor patients effectively today, which will help you become a better clinician tomorrow.

With soft-sealing eartips, adjustable headsets and a non-chill rim, these stethoscopes are comfortable for you and your patients. Whether you're examining adult or pediatric patients, the dual-sided chestpiece (with open or closed bell) features tunable diaphragms, so you can hear different frequency sounds without repositioning it. Choose the model right for you, each made in the USA without natural rubber latex.



3M™ Littmann® Cardiology IV™ Stethoscope

Helps isolate subtle sounds so you can hear changes in your patients in the most critical and challenging environments.



3M™ Littmann® Classic III™ Stethoscope

Compact and sensitive, perfect for assessing and monitoring in non-critical environments.

3M™ Littmann® Classic III™ Stethoscope



3M™ Littmann® Cardiology IV™ Stethoscope



Standard stainless steel



5620 Black 5621 Gray 5622 Navy 5623 Caribbean 5630 Ceil Blue



5631 Rose Pink 5633 Pearl Pink 5627 Burgundy 5831 Plum 5832 Lavender



5835 Turquoise 5839 Lemon-Lime 5648 Raspberry

Rainbow/brass finish

Rainbow/black

Copper finish



5806 Raspberry 5807 Caribbean 5870 Black 5809 Chocolate

Mirror/smoke



5959 Ceil Blue 5863 Navy 5960 Plum 5962 Pearl Pink

Champagne/smoke

Black/black



5861 Black 5864 Burgundy 5803 Black 5868 Burgundy

Smoke/smoke



5811 Black 5872 Turquoise 5873 Gray 5875 Lime Green

Product Features	
Acoustic Performance	7 (1-10 scale)
Chestpiece Design/Technology	Machined stainless steel, two-sided
Headset Material	Lightweight aluminum alloy
Length/Total Weight	69 cm/150 grams
Warranty	5 years
Diaphragm Details/Diameter	Two, single-piece tunable diaphragms 4,3 cm, 3,4 cm

Colors shown are approximate.

Explore Littmann stethoscopes at Littmann.co.uk



3M United Kingdom PLC
Charnwood Campus | 10 Bakewell Road
Loughbroug, Leicestershir, LE11 5RB

Standard stainless steel



6152 Black 6154 Navy 6155 Hunter 6156 Plum 6158 Raspberry



6184 Burgundy

Black/black



6163 Black 6168 Navy 6200 Black 6201 Black 6203 Black

Mirror

Rainbow/black



6177 Black 6170 Burgundy 6165 Black 6205 Plum

High-polish smoke/smoke

High-polish rainbow/smoke

High-polish rainbow/black



6238 Gray 6240 Black 6241 Raspberry 6242 Navy

High polish smoke/black



6202 Navy 6204 Black 6232 Black 6234 Caribbean

Champagne/black

Champagne/smoke

Champagne/champagne



6176 Burgundy 6190 Caribbean 6179 Black 6206 Hunter Green

Product Features	
Acoustic Performance	9 (1-10 scale)
Chestpiece Design/Technology	Machined stainless steel, two-sided
Headset Material	Lightweight aluminum alloy
Length/Total Weight	69cm/177 grams
Warranty	7 years
Diaphragm Details/Diameter	Two, single-piece tunable diaphragms 4,3 cm, 3,4 cm

3M, Littmann, the L Littmann logo, Classic III and Cardiology IV are trademarks of 3M. Please recycle. Printed in U.K. © 3M 2019. All rights reserved. 70-2011-7828-5. Rev. 06/2020